

Pendulum Lab using a Photogate Timer (instead of using a stopwatch)

NOTE: For *Pendulum Activity "Part 2: Data Collection and Analysis"*, the instruction sheet assumes that a stopwatch will be used to determine the period. Instead, we'll use a photogate (much easier!!!). So, you have to revise the data table structure, as shown below:

- a) *Period "T" (s) vs Mass "M" (kg)*, while keep length and displacement angle the same in every trial (L and θ are controlled)

Period "T" (s) vs Mass "M" (kg) of a Simple Pendulum

(Constant length, Displacement Angle less than 15 degrees)

1	2	3	4	5	6
	Trial 1 Period	Trial 2	Trial 3	(Average of the 3 trials)	average of 3 trials
Mass "M" (kg)	Time for 10 cycles "t"(s)	Time for 10 cycles "t"(s)	Time for 10 cycles "t"(s)	Time for 10 cycles "t"_{ave}(s)	Period "T"(s) = time for one cycle = t_{ave}/10 [column 5 divided by 10]
0.050					
0.100					
0.150					
0.200					
0.250					
0.300					
0.400					
0.500					
0.550					
0.600					

- b) *Period "T" (s) vs Length "L" (m)*, while keep mass and displacement angle the same in every trial (M and θ are controlled; θ should be kept small, at 15° or less)

Period "T" (s) vs Length "L" (m) of a Simple Pendulum

with a mass of 0.200 kg, and Displacement Angle less than 15 degrees

1	2	3	4	5	6	7	8
(ideal – try to get the string close to this length)	(measured actual value)	Trial 1	Trial 2	Trial 3	(Average of the 3 measured trials)	(Average of the 3 measured values)	(calculated from measured value)
		Period ←————→			Period		
Length "L" of pendulum (m) [ideal]	Length "L" of pendulum (m) (to the centre of mass of the pendulum)	Time for 10 cycles "t"(s)	Time for 10 cycles "t"(s)	Time for 10 cycles "t"_{ave}(s)	Time for 10 cycles "t"_{ave}(s)	Period, "T"(s) time for one cycle t_{ave}/10 [column 6 divided by 10]	Period squared "T" ² (s ²) [column 6 squared]
1.00							
0.90							
0.80							
0.70							
0.60							
0.50							
0.40							
0.30							
0.20							
0.10							

The Photogate timer

1. Use “Gate A” only (do not attach Gate B)
2. Set up the system as shown in this photo:



3. Set Mode to **[Period]**
4. Note that the pendulum has to swing through 3 cycles before the timer displays data.
5. As the pendulum swings for several cycles, watch the display and pay attention to the fluctuations in the data. *To what precision do you feel confident in the measurement?* (i.e. how many decimal places seem reasonably stable and constant?). The digits (numbers) that you are confident are stable in your measurement are called “**significant digits**” (or **significant figures = “sig figs”**). For example:
 - a. Pretend that as the pendulum swings through 5 cycles, you see the following readings for period: 0.6534 s; 0.6512 s; 0.6549 s; 0.6469 s; 0.6501 s.
 - b. You can see that the first 2 digits of the reading are stable and consistent (after rounding), but the last 2 digits fluctuate. In this case, you can feel confident that the period is $T = 0.65$ s, but the 3rd and 4th decimal places are uncertain. On the data table, you should include only the “certain” digits. You would write **0.65 s** in the column for period. This measurement is precise to the nearest 100th of a second, with 2 significant digits.
6. In the data table column for “Period” (T), write only the digits that you are confident are stable and correct (i.e. include only a reasonable and appropriate number of **significant digits**).